

# Discrimination of Proximal Hip Fracture by Quantitative Ultrasound Measurement at the Radius

M. Weiss, A. Ben-Shlomo, P. Hagag, and S. Ish-Shalom  
Osteoporosis Int. (2000) 11:411-416

Omnisense demonstrates its ability to assess hip fracture risk of subjects by speed of sound measurements at the radius.

## Introduction

Osteoporosis is a disease that culminates with fragility fractures and imposes a major burden on health care costs. In dealing with this worldwide condition, it is prudent to use a reliable, non-expensive, portable diagnostic device that does not use ionizing radiation and is capable of measuring bone properties at several sites. Sunlight Omnisense 7000S, a quantitative ultrasound device which measures the speed of sound (SOS) along the bone at multiple skeletal sites, provides such a solution.

## Study Design

The study was designed to evaluate the ability of the Omnisense to assess fracture risk.


Three groups of women were recruited for this cross-sectional study:

- Elderly women with osteoporotic proximal hip fractures
- Elderly women without any evidence of osteoporotic fractures
- Young, healthy women

Speed of Sound measurements were performed at the distal 1/3 of the radius.

## Results

SOS was found to be significantly lower in the fractured group than in the non-fractured group. Hip fracture subjects had a mean speed of sound (SOS) of  $3861 \pm 149$  m/s, while non-fracture subjects had a mean SOS of  $3966 \pm 145$  m/s. Thus, there was a difference of 105 m/s between the two groups (t-test  $p < 0.0001$ ). A low SOS is an indication of unhealthy bone and increased fracture risk (ROC curves indicate an AUC of 0.79 (95% CI, 0.73 - 0.86)). Assessment of fracture risk can be presented by odds ratio (OR). In this study, OR was 2.16 (95% CI, 1.46 - 3.19). This analysis shows that for every 100 m/s decrease in SOS the odds of fracture increase by about 50% and that for every decrease of 162 m/s in SOS the odds of fracture double.



## Conclusion ►►

SOS, as measured by the Omnisense, can be considered an important factor in aiding the physician when diagnosing a patient for osteoporosis and determining the patient's risk of fracture.