

SuperSonic Imagine's Innovative Solutions for Clinical Liver Disease

Global Product Marketing & Education
October 2016




Liver Disease is a Growing Global Problem

The Facts:

- Liver fibrosis is a GLOBAL problem impacting 300-700 million people (5-10% of the world's population).
- The most common causes are:


Hepatitis B virus (HBV), Hepatitis C (HBV), Alcohol and Non-alcoholic Fatty Liver Disease (NAFLD) including Non-alcoholic Steatohepatitis (NASH).



Reference: The global impact of hepatic fibrosis and end-stage liver disease.
Lim YS1, Kim WR. Clin Liver Dis. 2008 Nov;12(4):733-46, vii. doi: 10.1016/j.cld.2008.07.007.


Image: The global burden of liver disease: a challenge for methods and for public health
Peter Byass, BMC Medicine 2014 12:159.

Global Prevalence of HBV



240 million people worldwide


Global Prevalence of HCV




180 million people worldwide

Prevalence of HCV per Country

| Country | Income classification | Most prevalent genotypes | Anti-HCV (%) | No. infected |
|------------------------------|-----------------------|--------------------------|--------------|--------------|
| China | Upper-middle | 1,2,6 | 2.2 | 29,791,212 |
| India | Lower-middle | 1,3 | 1.5 | 18,216,960 |
| Egypt | Lower-middle | 4 | 14 | 11,826,360 |
| Indonesia | Lower-middle | 1,2 | 3.9 | 9,436,986 |
| Pakistan | Lower-middle | 3 | 5.9 | 9,422,403 |
| Russia | Upper-middle | 1,3 | 4.1 | 5,796,498 |
| USA | High | 1,2,3 | 1.8 | 5,367,834 |
| Democratic Republic of Congo | Low | 4 | 6.4 | 4,010,240 |
| Nigeria | Lower-middle | 1,2 | | |
| Japan | High | 1,2 | | |
| Cameroon | Lower-middle | 1,2,4 | | |
| Brazil | Upper-middle | 1,3 | | |
| Uganda | Low | 1,4 | | |
| Philippines | Lower-middle | 1 | | |
| Italy | High | 1,2 | | |
| Ukraine | Lower-middle | 1 | | |
| Uzbekistan | Lower-middle | 1,3 | | |
| Turkey | Upper-middle | 1 | | |
| Ethiopia | Low | 1,2,4 | | |
| Thailand | Upper-middle | 1,3,6 | | |
| World's Population | | | | |





The Clinical Challenge: What do clinicians want?



A faster, easier, safer and non-invasive alternative to liver biopsy.

A method which allows for repeated serial evaluations. This is critical for monitoring new anti-viral HCV drugs and NAFLD-NASH.



A validated tool to screen for liver fibrosis/cirrhosis, AND the ability to stage fibrosis with a score.

A method which is low cost for the patient, AND provides a good reimbursement and return on investment.

- Current costs: Up to \$300,000 over a patients lifetime



There is global consensus that a non-invasive, accurate, validated, and economical test for early detection and staging is needed.

Our Solution:
The Aixplorer[®] ultrasound system
for Liver Fibrosis



Aixplorer's Key Pillars




“One Probe Solution”
Unique ShearWave Technology



“60 Second Exam”
Faster Workflow



“Clinically Proven”
Superior Performance and Accuracy



“Higher Reimbursement”
4 times Transient Elastography

Now there is a Solution

Aixplorer® with ShearWave™ Elastography



Our Solution:

- **SuperSonic Imagine's** unique technology can provide liver fibrosis **screening, staging and monitoring.**
- The 60 second ShearWave™ Elastography exam – **fast, easy, clinically proven** and **cost effective:**
 - Simple 1 image acquisition and 3 measurement workflow
 - Quick Q-Box™ Analysis Tools with Stability Index to assure quality
 - Clear and concise Liver Report package



One Probe Solution

XC6-1 Single Crystal Curved Array Transducer




- “One Probe Solution” for All Patient Types:

- Reliable – One probe lasts a lifetime. No need for annual recalibration.
- Excellent SWE™ penetration down to 10 cm in obese patients. Maximum imaging to depths beyond 30 cm!
- No failures due to ascites.
- Single crystal technology for maximum sensitivity and image clarity.
- Lightweight, ergonomic design.




60 Second Exam

Non-Invasive & Quantitative Liver Fibrosis Assessment



- Q-Box™ quantitative tools allow rapid quantitative measurement of tissue stiffness.
- Stability Index for accurate and confident quantitative results.

- Large color-coded ShearWave™ area for easy visualization of the heterogeneity and extent of liver fibrosis.
- Simple blue to red color scale for easy qualitative interpretation of stiffness.




60 Second Exam

4 Easy Steps

1. Position the patient in a supine or slight left oblique position, hand above head, and activate the XC6-1 probe.
2. Firmly press the probe to enlarge the intercostal space. When a clear B-mode image of the liver is seen up to 8 cm of depth, activate SWE™ and ask the patient to pause breathing.
3. Acquire 5-10 seconds of real-time SWE images. Press Freeze.
4. Choose an image with clear parenchyma free of vessels. Press Q-Box™. Position the ROI in an area where the Stability Index (SI) reads 90% or higher. Press Save.

Repeat this procedure 3 times to acquire 3 valid, independent SWE™ images of the same scanning view.



The Aixplorer[®]: Clinically Proven for Liver Fibrosis



Clinically Proven

Liver Fibrosis Assessment with SWE™



Liver stiffness values measured with SWE in Hepatitis C patients


| METAVIR Scores | Stiffness values (IQR) | Fibrosis assessment (METAVIR) | Cut-off values |
|-----------------------|-------------------------------|--------------------------------------|-----------------------|
| F0-F1 | 5.1-6.8 | | |
| F2 | 7.2-8.3 | F \geq 2 | 7.1 |
| F3 | 9.2-10.1 | F \geq 3 | 8.7 |
| F4 | 12.8-18.8 | F=4 | 10.4 |

Values taken from “Accuracy of real-time shear wave elastography for assessing liver fibrosis in chronic hepatitis C: a pilot study. Ferraioli G, Tinelli C, Dal Bello B, Zicchetti M, Filice G, Filice C; Liver Fibrosis Study Group. Hepatology. 2012 Dec;56(6):2125-33.”


Clinically Proven Liver Fibrosis Assessment with SWE™



HCV Patients



NAFLD Patients



Bavu E et al. Ultrasound Med Biol. 2011 Sep;37(9):1361-73.

Ferraioli G et al. Hepatology. 2012 Dec;56(6):2125-33.

Yoneda M et al. Clin Gastroenterol Hepatol. 2015 Aug;13(8):1502-9.

(% Patients with BMI > 25 kg/m²)

Herrmann et al. J Hepatol 2015 Apr;62:S187–S212. (£ Oral presentation at EASL 2015)

Herrmann et al. J Hepatol 2015 Apr;62:S187–S212. (£ Oral presentation at EASL 2015)

Clinically Proven

Liver Fibrosis Assessment with SWE™




- Correlation coefficient between liver stiffness measured non-invasively and liver fibrosis scoring

| | SWE™ | TE | ARFI VTTQ |
|-------------------------------|-------------|-------------|-------------|
| Ferraioli¹ | 0.83 | 0.74 | |
| Cassinotto² | 0.79 | 0.70 | 0.64 |
| Gerber³ | 0.71 | 0.73 | 0.75 |
| Samir⁴ | 0.41 | | |
| Dhyani⁵ | 0.58 | | |
| Franchi⁶ | 0.83 | | |

1 Ferraioli G et al. Hepatology. 2012 Dec;56(6):2125-33.
 2 Cassinotto C et al. J Hepatol. 2014 Sep;61(3):550-7.
 3 Gerber L et al. Ultrasound Med Biol. 2015 Sep;41(9):2350-9.
 4 Samir A et al. Radiology. 2015 Mar;274(3):888-96.
 5 Dhyani M et al. J Med Imaging Radiat Oncol. 2015 Dec;59(6):687-94.
 6 Franchi-Abella S et al. Radiology. 2016 Feb;278(2):554-62.

SSI SWE: No overlap between IQRs



Clinically Proven SWE™ Multicenter Retrospective Study



First Results of the Multicenter Retrospective Study on Liver Fibrosis Assessment with SWE™

15 International Sites / 1340 patients

ShearWave™ Elastography delivers **excellent** results in managing chronic liver disease patients

In HBV patients, SWE™ delivers **increased sensitivity** in the assessment of significant fibrosis and cirrhosis ($\leq F1$ vs. $\geq F2$)

In HCV patients, SWE delivers **excellent sensitivity** in the evaluation of cirrhosis ($\geq F4$)

The Results Are IN!

RS-3186
Non-Invasive markers of liver fibrosis

2D-SHEAR WAVE ELASTOGRAPHY IS EQUIVALENT OR SUPERIOR TO TRANSIENT ELASTOGRAPHY FOR LIVER FIBROSIS ASSESSMENT:

RESULTS FROM AN INDIVIDUAL PATIENT DATA BASED META-ANALYSIS

Eva Herrmann¹, Victor de Lédinghen², Christophe Cassinotto³, Winnie C.-W. Chu⁴, Vivian Y.-F. Leung⁵, Giovanna Ferraioli⁶, Carlo Filice⁶, Laurent Castera⁷, Valérie Vilgrain⁸, Maxime Ronot⁸, Jérôme Dumortier⁹, Aymeric Guibal¹⁰, Stanislas Pol¹¹, Jonel Trebicka¹², Christian Strassburg¹², Ronqin Zhena¹³, Jian Zhena¹³, Sven Vonghia¹⁴, Emanuel K. Manesis¹⁵, Pavlos Zoumpoulis¹⁶, Ioan Sporea¹⁷

ShearWave™ Elastography delivers **excellent diagnostic** assessment of all fibrosis stages

Department of Ultrasound Imaging, 3rd Hospital of Sun Yat-Sen University, Guangzhou, China, ¹⁴Division of Gastroenterology and Hepatology, Antwerp University Hospital, Edegem, Belgium, ¹⁵Department of Internal Medicine, ¹⁶Department of Radiology, Athens University School of Medicine, Athens, Greece, ¹⁷Department of Gastroenterology and Hepatology, Victor Babes University of Medicine and Pharmacy, Timisoara, Romania, ¹⁸Research Unit for Gastroenterology and Hepatology, Odense University Hospital, Odense, Denmark, ¹⁹Department of Internal Medicine, Goethe University Hospital Frankfurt, Frankfurt, Germany

Corresponding author's email: herrmann@med.uni-frankfurt.de

Do you want to apply for a Young Investigator Bursary?: No

Background and Aims: 2D shear wave elastography (2D-SWE) based on supersonic shear imaging (SSI) has proven to be efficient for liver fibrosis assessment in several small to moderate size trials. We aimed at running a larger scale analysis.

Methods: Data from 13 sites were collected retrospectively from 1340 patients. The database was cleared by the French National Commission for Data Protection. The data were analysed centrally.

Results: The study included 1340 patients with liver disease of various etiologies. The database was cleared by the French National Commission for Data Protection. The data were analysed centrally.

Overall, 19.3% had significant fibrosis, 11.0% had severe fibrosis, and 3.7% had cirrhosis. The database was cleared by the French National Commission for Data Protection. The data were analysed centrally.

Optimal cut-offs were 7.1 kPa for diagnosing significant fibrosis in all patients (75.7% correctly classified), 13.5 kPa for diagnosing cirrhosis in HCV and NAFLD patients, and 11.5 kPa for diagnosing cirrhosis in HBV patients (87% correctly classified).

Differences in AUROC were borderline significant for diagnosing significant fibrosis (95%>CI for AUROC-2D-SWE minus AUROC-TE: [0.0004, 0.055], p=0.047) and AUROC was significantly higher for 2D-SWE when diagnosing cirrhosis (95%>CI for AUROC-2D-SWE minus AUROC-TE: [0.006, 0.036], p=0.0058). 2D-SWE was superior for diagnosing

Clinically Proven Liver Fibrosis Assessment with SWE™ vs. FibroScan®



Quantitative Elastography of Liver Fibrosis and Spleen Stiffness in Chronic Hepatitis B Carriers: Comparison of Shear-Wave Elastography and Transient Elastography with Liver Biopsy Correlation¹

Purpose: To document utility of shear-wave (SW) elastography assessing liver fibrosis in chronic hepatitis B and to compare SW elastography with transient elastography (TE) and liver biopsy.

Conclusion: SW elastography provides more accurate liver elasticity with liver fibrosis stage than TE and liver biopsy, especially in intermediate and advanced stages.

Vivian Yee-fong Leung, PhD

ShearWave™ Elastography has Higher Accuracy and Technical Success than the FibroScan® in Hepatitis B patients

¹From the Institute of Digestive Disease (J.S., V.W.W., G.L.W., A.M.C., S.H.C., H.L.C., W.C.C.), Department of Imaging and Interventional Radiology (V.Y.L., J.A., A.T.A., W.C.C.), Department of Medicine and Therapeutics (J.S., V.W.W., G.L.W., A.M.C., S.H.C., H.L.C.), and Department of Anatomical and Cellular Pathology (A.W.C., P.C.C.), Prince of Wales Hospital, The Chinese University of Hong Kong, Ngan Shing Street, Shatin, Hong Kong, SAR, China 852. Received January 16, 2013; revision requested February 22; revision received March 20; accepted May 13; final version accepted May 23. Address correspondence to W.C.C. (e-mail: winnie@med.cuhk.edu.hk).

© RSNA, 2013

910

radiology.rsna.org • Radiology: Volume 269: Number 3—December 2013

Radiology

CLINICAL—LIVER

Transient and 2-Dimensional Shear-Wave Elastography Provide Comparable Assessment of Alcoholic Liver Fibrosis and Cirrhosis

Maja Thiele,^{1,2,3} Sönke Detlefsen,^{3,4} Linda Sevelsted Møller,⁵ Bjørn Stæhr Madsen,^{1,2,3} Janne Fuglsang Hansen,^{3,6} Annette Dam Fialla,^{1,3,7} Jonel Trebicka,⁸ and Aleksander Krag^{1,3}

¹Department of Gastroenterology and Hepatology, Odense University Hospital, Odense, Denmark; ²OPEN Odense Liver Disease Collaborative Network, Odense University Hospital, Odense, Denmark; ³Institute of Clinical Research, University of Southern Denmark, Odense, Denmark; ⁴Department of Pathology, Odense University Hospital, Odense, Denmark; ⁵Department of Medicine, Odense University Hospital, Svendborg, Denmark; ⁶Department of Infectious Diseases, Odense University Hospital, Odense, Denmark; ⁷Department of Medicine, Hospital of Southwest Jutland, Esbjerg, Denmark; and ⁸Department of Medicine I, University of Bonn, Bonn, Germany

See editorial on page 29.

BACKGROUND & AIMS: Alcohol abuse causes half of all deaths from cirrhosis in the West, but few tools are available for noninvasive diagnosis of alcoholic liver disease. We evaluated 2 elastography techniques for diagnosis of alcoholic liver disease and cirrhosis: liver biopsy with Ishak score and collagen-proportionate area were used as reference. **METHODS:** We performed a prospective study of 199 consecutive patients with ongoing or prior alcohol abuse, but without known liver disease. One group of patients had a high pretest probability of cirrhosis because they were identified at hospital liver clinics (in Southern Denmark). The second, lower-risk group, was recruited from municipal alcohol rehabilitation centers and the Danish national public health portal. All subjects underwent same-day transient elastography (FibroScan). 2-dimensional shear wave elastography

ShearWave™ Elastography yields narrower IQR ranges and tighter probability curves in the risk prediction of alcoholic fibrosis and cirrhosis

at risk for liver fibrosis due to alcohol consumption, we found elastography to be an excellent tool for diagnosing liver fibrosis and for excluding (ruling out rather than ruling in) cirrhosis.

Keywords: Supersonic Shear Imaging; AUC; Noninvasive Methods; Diagnostic Test.

Abbreviations used in this paper: 2D-SWE, real-time 2-dimensional shear wave elastography; ALD, alcoholic liver disease; AUC, area under the receiver operating characteristics curve; CI, confidence interval; CPA, collagen proportionate area; GGT, gamma-glutamyltransferase; IQR, interquartile range; TE, transient elastography.

Most current article

© 2016 by the AGA Institute
0016-5085/\$36.00

<http://dx.doi.org/10.1053/j.gastro.2015.09.040>

HEPATOLOGY

Official Journal of the American Association for the Study of Liver Diseases



Accuracy of Real-Time Shear Wave Elastography for Assessing Liver Fibrosis in Chronic Hepatitis C: A Pilot Study

Anna Ferraioli,¹ Carmine Tinelli,² Barbara Dal Bello,³ Mabel Zicchetti,¹ Gaetano Filice,⁴ and Carlo Filice¹ on behalf of the Liver Fibrosis Study Group

Real-time shear wave elastography (SWE) is a novel, noninvasive method to assess liver fibrosis by measuring liver stiffness. This single-center study was conducted to assess the accuracy of SWE in patients with chronic hepatitis C (CHC), in comparison with transient elastography (TE), by using liver biopsy (LB) as the reference standard. Consecutive patients with CHC scheduled for LB by referring physicians were studied. One hundred twenty-one patients met inclusion criteria. On the same day, real-time SWE using the sound (US) system, Aixplorer (SuperSonic Imagine S.A., Aix-en-Provence, France),

ShearWave™ Elastography is More Accurate than the FibroScan® in Hepatitis C patients

00:000-000

Chronic hepatitis C (CHC), prognosis and management largely by the extent of fibrosis. Liver biopsy (LB) is still considered the gold standard for evaluation of liver fibrosis, even though it is painful, costly, and with limitations in accuracy. The accuracy of LB is influenced by many factors, as well as intra- and interobserver error. Given these limitations, LB is not ideal for repeated assessment of disease

progression. Following not only the progression, but also the regression of liver fibrosis over time could be of clinical significance, because research has demonstrated reduction in liver fibrosis with treatment, even in advanced stages.^{5,6}

These limitations of the LB have motivated research for noninvasive methods of measuring liver fibrosis. Transient elastography (TE) has emerged as the noninvasive test of reference and is entering clinical practice in Europe.⁷ TE is a noninvasive method that evaluates

Conclusion: Real-time SWE using the Aixplorer system was available for the study by SuperSonic Imagine S.A. (Aix-en-Provence, France).

1

Clinically Proven Liver Fibrosis Assessment with SWE™ vs. FibroScan®



Author's personal copy

Research Article

Non-invasive assessment of liver fibrosis with impulse elastography: Comparison of Supersonic Shear Imaging with ARFI and FibroScan®

Christophe Cassinotto^{1,2,*}, Bruno Lapuyade¹, Amaury Mouries¹, Jean-Baptiste Hiriart³, Julien Vergniol³, Delphine Gaye¹, Claire Castain¹, Brigitte Le Bail^{2,4}, Faiza Chermak³, Juliette Foucher³, François Laurent⁵, Michel Montaudon⁵, Victor De Ledinghen^{2,3}

¹Department of Diagnostic and Interventional Imaging, Hôpital Haut-Lévêque, Centre Hospitalier Universitaire de Bordeaux, 1 Avenue de Magellan, 33604 Pessac, France; ²INSERM U1053, Université Bordeaux, Bordeaux, France; ³Centre d'investigation de la fibrose hépatique, Hôpital Haut-Lévêque, Centre Hospitalier Universitaire de Bordeaux, 1 Avenue de Magellan, 33604 Pessac, France; ⁴Department of Pathology, Hôpital Pellegrin, Centre Hospitalier Universitaire de Bordeaux, Place Amélie Raba-Léon, 33000 Bordeaux, France; ⁵Department of Diagnostic and Interventional Imaging, Hôpital Haut-Lévêque, CHU and University of Bordeaux, 1 Avenue de Magellan, 33604 Pessac, France

Background & Aims: Non-invasive assessment of liver fibrosis by elastography is a rapidly developing field with frequent technological innovations. The aim of this study was to assess the diagnostic performances of Supersonic Shear Imaging (SSI) for the diagnosis of liver fibrosis in chronic liver disease.

Methods: A total of 349 consecutive patients with chronic liver diseases who underwent liver biopsy from November 2011 to October 2013 were prospectively enrolled. For each patient, liver stiffness was assessed by SSI, ARFI, FibroScan® (M probe for patients with BMI <30 kg/m², and XL probe for patients with

than ARFI for the diagnosis of significant fibrosis (>F2) (p = 0.0003). No significant difference was observed for the diagnosis of mild fibrosis and cirrhosis.

Conclusions: SSI is an efficient method for the assessment of liver fibrosis in chronic liver diseases, comparing favourably to FibroScan® and ARFI.

© 2014 European Association for the Study of the Liver. Published by Elsevier B.V. All rights reserved.

Clinical Gastroenterology and Hepatology 2015;13:1502-1509

Supersonic Shear Imaging and Transient Elastography With the XL Probe Accurately Detect Fibrosis in Overweight or Obese Patients With Chronic Liver Disease

Masato Yoneda, Emmanuel Thomas, Seth N. Sclair, Tiffannia T. Grant, and Eugene R. S. Schiff

Schiff Center for Liver Diseases, University of Miami Miller School of Medicine, Miami, Florida

BACKGROUND & AIMS: Assessment of the severity of liver fibrosis is an important step in evaluating patients with chronic liver disease and determining their prognosis. We compared liver stiffness measurements (LSMs) made by supersonic shear imaging (SSI) with those of transient elastography (TE)-XL for their ability to determine the degree of liver fibrosis in overweight or obese patients with chronic liver disease.

ShearWave™ Elastography is Comparable to Biopsy and Superior to FibroScan® in Obese Patients

the 102 biopsy-proven patients with chronic HCV infection.

CONCLUSIONS: SSI and the TE-XL probe each accurately quantify liver fibrosis in overweight or obese patients with chronic liver disease, including those with HCV infection, when compared with data obtained from histologic analysis. SSI data obtained from the right lobe and the TE-XL probe can be used to evaluate fibrosis with similar accuracy.

Keywords: Viral Hepatitis; Cirrhosis; Fibrosis Stage; Obesity.

See editorial on page 1510.

Staging of liver fibrosis is essential in determining the prognosis and optimal treatment for patients with chronic liver disease and also to guide surveillance for the development of hepatocellular carcinoma.¹ Liver biopsy is recommended as the reference standard method for the diagnosis and staging of fibrosis in chronic liver disease.² This procedure, however, is costly, and is time consuming both for providers and patients.³ In addition, despite being the gold standard test for assessing liver fibrosis, liver biopsy is limited

further by sampling error and intra-observer/inter-observer variability.^{4,5} Therefore, there is a need for rapid, quantitative, and noninvasive methods for the assessment of liver fibrosis.

Abbreviations used in this paper: AAR, aspartate aminotransferase ratio; APRF, aspartate aminotransferase to platelet ratio index; AST, aspartate aminotransferase; ALT, alanine aminotransferase; BMI, body mass index; F4, fibrosis stage 4; FIB-4, Fibrosis-4 score; HCVR, hepatic compliance ratio; LSM, liver stiffness measurement; SSI, supersonic shear imaging from the left lobe; SSI-R, supersonic shear imaging from the right lobe; TE, transient elastography.

© 2015 by the AGA Institute
1542-3565/\$36.00
<http://dx.doi.org/10.1016/j.cgh.2015.03.016>

ShearWave™ Elastography Clinically Outperforms the FibroScan® in Chronic Liver Disease Patients

the measurement of the velocity of a local shear wave through soft tissues. As ARFI and conversely to FibroScan®, this method is built on an ultrasound device (Aixplorer, Supersonic Imagine, Aix-en-Provence, France), and requires no external vibrator to generate the shear wave.

Journal of Hepatology 2014 vol. 61 | 550-557

ELSEVIER

Prospective Comparison of Spleen and Liver Stiffness by Using Shear-Wave and Transient Elastography for Detection of Portal Hypertension in Cirrhosis¹

Purpose: To prospectively compare the technical success rate and accuracy of shear-wave elastography (SWE) and transient elastography (TE) for the detection of clinically significant portal hypertension (PH) in patients with advanced cirrhosis who are undergoing hepatic vein pressure gradient (HVPG) measurements.

Materials and Methods: The institutional ethics committee approved the study, and written informed consent was obtained. Seventy-nine consecutive patients with cirrhosis who were undergoing SWE and TE at the time of HVPG measurement were studied. The technical success rate of SWE and TE was compared with the diagnostic value of liver stiffness (LS) and spleen stiffness (SS) measurements and composite scores (LS spleen-diameter-to-platelet-ratio score [LSPS] and PH risk score) by using SWE and TE to detect clinically significant PH (HVPG ≥ 10 mm Hg) and esophageal varices at high risk of bleeding. Areas under the receiver operating characteristic curve and the DeLong test were

Laurent Castéra, MD, PhD
Pierre-Emmanuel Rautou, MD, PhD
Maxime Ronot, MD, PhD
Simon Lambert, PhD
Marco Dioguardi Burgio, MD
Claire Francoz, MD, PhD
Urbelle Plessier, MD
François Durand, MD
Dominique Valla, MD
Didier Lebrec, MD
Alérie Vilgrain, MD, PhD
Laurent Castéra, MD

ShearWave™ Elastography is more Technically Successful and has Better Performance than the FibroScan® in the Detection of Portal Hypertension in Cirrhosis

was significantly better than that of FibroScan® (97% vs 44%, P < .001). The technical success rate of SWE was significantly better than that of FibroScan® (97% vs 44%, P < .001). The diagnostic value of SWE was significantly better than that of FibroScan® (97% vs 44%, P < .001).

© RSNA, 2014

Radiology: Volume 000: Number 0—■■■ 2015 • radiology.rsna.org

Clinically Proven

What the experts are saying



“I can reduce my biopsy numbers by 90% during initial fibrosis staging when I use the Aixplorer’s® real-time, quantitative ShearWave™ Elastography from SuperSonic Imagine”

“This is especially helpful for my patients who are eligible to receive the new antiviral treatments for HCV. This quick, non-invasive exam improves the overall patient experience.”

Dr. Ravi Ravinuthala
Ohio Gastroenterology & Liver Institute
Cincinnati, Ohio

Dr. James Trotter
Medical Director of Liver Transplantation
Baylor University Medical Center
Dallas, Texas



“The advent of ShearWave™ Elastography, has ushered in an era of fewer liver biopsies for the Hepatitis C patient, while allowing the clinician to evaluate, monitor and effectively treat these patients without the risks of an invasive biopsy.”


The Clinical Challenge: Liver Steatosis



The Clinical Challenge

Assessment of Liver Steatosis

- What is Steatosis?
 - Steatosis, or fatty liver, describes the accumulation of excess fat in the liver, usually above 5%.
 - Steatosis is the hallmark of most prevalent liver diseases, including alcoholic and nonalcoholic fatty liver disease (NAFLD).



- Steatosis / NAFLD > 5% may progress to nonalcoholic steatohepatitis (NASH), fibrosis, cirrhosis and liver cancer.
- Steatosis also reduces the viral response in hepatitis C treatment, complicates liver surgery, and raises overall risk to cardiovascular events.

The Clinical Challenge


Obesity and Steatosis



The Facts:

- NAFLD and NASH, and their progression was relatively unknown in etiology until 10 years ago.
- Today, more than 25% of the world's population is at risk with NAFLD, 3% with NASH.
- NAFLD/NASH has very high prevalence in the Americas, Asia-Pacific, the Middle East and Europe.
- Approximately 10%-25% of patients with NAFLD will develop NASH. 5%-8% of those will develop liver cirrhosis within 5 years. Furthermore, 12.8% of patients with liver cirrhosis will develop hepatocellular carcinoma (HCC) within 3 years.
- NAFLD/NASH is a silent killer, growing at an epidemic rate in proportion to obesity. From 1980-2013 obesity increased by 10% in adults and by 47% in children.

| Region | Population studied | Prevalence of NAFLD in these populations (%) |
|-------------------|----------------------|--|
| USA | Pediatric population | 13–14 |
| | General population | 27–34 |
| | Morbid obesity | 75–92 |
| | European-Americans | 33 |
| | Hispanic-Americans | 45 |
| | African-Americans | 24 |
| Europe | Pediatric population | 2.6–10 |
| | General population | 20–30 |
| Western countries | General population | 20–40 |
| | Obesity or diabetes | 75 |
| | Morbid obesity | 90–95 |
| Worldwide | Obese population | 40–90 |
| Middle East | General population | 20–30 |
| Far East | General population | 15 |
| Pakistan | General population | 18 |





Non-alcoholic fatty liver disease and obesity: Biochemical, metabolic and clinical presentations. Milić S. World Journal of Gastroenterology : WJG. 2014;20(28):9330-9337. doi:10.3748/wjg.v20.i28.9330.

USA Predilection for Liver Fibrosis

USA Statistics:

- 700,000-1.4 million persons are estimated to be infected with the HBV virus
- 3.2 million persons are chronically HCV-infected
- 27-34% of all Americans are estimated to have NAFLD



| Region | Population studied | Prevalence of NAFLD in these populations (%) |
|--------|----------------------|--|
| USA | Pediatric population | 13-14 |
| | General population | 27-34 |
| | Morbid obesity | 75-92 |
| | European-Americans | 33 |
| | Hispanic-Americans | 45 |
| | African-Americans | 24 |

USA Data: <http://www.cdc.gov/hepatitis/statistics/2013surveillance/commentary.htm>

The Clinical Challenge

Obesity and Steatosis



“Left unchecked, obesity will make the current generation of children the first in human history to have a life span shorter than that of their parents.”

David Satcher, MD, PhD
US Surgeon General (1998-2002)



A Global Effort to Address the Epidemic of Fatty Liver Disease
William F. Balistreri. Medscape. September 24, 2014

The Clinical Challenge

Assessment of Liver Steatosis



Low cost.



Non-invasive.
Safe.
Well tolerated.



Early detection
of steatosis at
the 5% level.



Serial assessment
during treatment
and beyond.



Globally
available.

A non-invasive, accurate, validated, and economical test to detect, assess and monitor steatosis is needed.

Our Solution:
The Aixplorer[®] ultrasound system
for Liver Steatosis



60 Second Exam

Steatosis Assessment with B-mode Ratio Tool



- Hepatorenal Ratio for Steatosis Assessment: B-Mode Ratio Tool
 - Quantifies brightness of the liver with renal cortex.
 - Significantly improves the assessment of fatty liver compared to traditional “guess-timation”.
 - Fast and easy to use. Requires only B-mode.
- Clinically Proven
 - Ultrasound Hepatic-Renal Ratio was found to be an excellent predictor of liver fat content.
 - Able to discriminate between levels of steatosis with high sensitivity and specificity.^[1]
 - The addition of ultrasound attenuation parameter (e.g. CAP) only improved the accuracy by 1.8%.^[2]
 - Can be done on DICOM images.^[3]
- No additional cost
 - Compare to FibroScan CAP tool at +\$22k.

Gastrointestinal Imaging
Original Research

Diagnostic Value of a Computerized Hepatorenal Index for Sonographic Quantification of Liver Steatosis

Muriel Webb¹
Hanny Yeshua^{1,2}
Shira Zelber-Sagi¹
Erwin Santo¹
Eli Brazowski³
Zamir Halpern¹
Ran Oren¹

OBJECTIVE. Quantification of liver steatosis is clinically relevant in various liver diseases but cannot be done by conventional sonography, which only provides a qualitative assessment with significant observer variability. The aim of this study was to assess sonography as an objective tool for the quantification of liver steatosis.

MATERIALS AND METHODS. Files of 111 patients with chronic liver disease who were referred for sonographically guided liver biopsy were collected. A hepatorenal sonographic index was calculated on the basis of the ratio between the echogenicity of the liver and that of the right kidney cortex using histogram echo intensity. Liver steatosis was graded

| +B-mode Ratio | |
|---------------|--------|
| B Ratio | 1.09 |
| Mean | 85.8 |
| Diam | 6.0 mm |
| Depth | 6.6 cm |
| Mean | 78.9 |
| Diam | 6.0 mm |
| Depth | 6.6 cm |

Normal < 1.5

AJR:192, April 2009 909

[1] Diagnostic value of a computerized hepatorenal index for sonographic quantification of liver steatosis. Webb M et al. AJR Am J Roentgenol. 2009 Apr;192(4):909-14.

[2] Standardized Ultrasound Hepatic/Renal Ratio and Hepatic attenuation Rate to Quantify Liver Fat Content: an Improvement Method. Xia, et. al. Obesity (2011) 20, 444–452.

[3] Improved method for calculating hepatic steatosis using the hepatorenal index. Shiralkar K et al. J Ultrasound Med. 2015 Jun;34(6):1051-9.

The Aixplorer[®]: Clinically Proven for Liver Steatosis



Clinically Proven Ultrasound Meta-Analysis for the Detection of Steatosis



Diagnostic Accuracy and Reliability of Ultrasonography for the Detection of Fatty Liver: A Meta-Analysis

49 Publications / 4720 patients

Ultrasonography allows for
reliable and accurate detection
of moderate-severe fatty liver,
compared to histology.

The Hepato-Renal brightness (HRI)
performed best
with a sensitivity of 98% and specificity of 93%.

Other measures of ultrasound brightness
e.g. hepatic or portal vessel walls
significantly outperformed
ultrasound beam attenuation (CAP),
sensitivity 59%, specificity 95%.

HEPATOLOGY, Vol. 54, No. 3, 2011

Diagnostic Accuracy and Reliability of Ultrasonography for the Detection of Fatty Liver: A Meta-Analysis

Ruben Hernaez,^{1,2,3*} Mariana Lazo,^{1*} Susanne Bonekamp,⁴ Ihab Kamel,⁴ Frederick L. Brancati,^{1,3,5}
Eliseo Guallar,^{3,5,6} and Jeanne M. Clark^{1,3,5}

Ultrasonography is a widely accessible imaging technique for the detection of fatty liver, but the reported accuracy and reliability have been inconsistent across studies. We aimed to perform a systematic review and meta-analysis of the diagnostic accuracy and reliability of ultrasonography for the detection of fatty liver. We used MEDLINE and Embase from October 1967 to March 2010. Studies that provided cross-tabulations of ultrasonography versus histology or standard imaging techniques, or that provided reliability data for ultrasonography, were included. Study variables were independently abstracted by three reviewers and double checked by one reviewer. Forty-nine (4720 participants) studies were included for the meta-analysis of diagnostic accuracy. The overall sensitivity, specificity, positive likelihood ratio, and negative likelihood ratio of ultrasound for the detection of moderate-severe fatty liver, compared to histology (gold standard), were 84.8% (95% confidence interval: 79.5-88.9), 93.6% (87.2-97.0), 13.3 (6.4-27.6), and 0.16 (0.12-0.22), respectively. The area under the summary receiving operating characteristics curve was 0.93 (0.91-0.95). Reliability of ultrasound for the detection of fatty liver showed kappa statistics ranging from 0.54 to 0.92 for intrarater reliability and from 0.44 to 1.00 for interrater reliability. Sensitivity and specificity of ultrasound was similar to that of other imaging techniques (i.e., computed tomography or magnetic resonance imaging). Statistical heterogeneity was present even after stratification for multiple clinically relevant characteristics. **Conclusion:** Ultrasonography allows for reliable and accurate detection of moderate-severe fatty liver, compared to histology. Because of its low cost, safety, and accessibility, ultrasound is likely the imaging technique of choice for screening for fatty liver in clinical and population settings. (HEPATOLOGY 2011;54:1082-1090)

Fatty liver is the accumulation of fat (i.e., macrovesicular steatosis) within the hepatic parenchyma. Nonalcoholic fatty liver disease (NAFLD), the presence of fat infiltration in the liver in the absence of excessive alcohol consumption and other causes of liver disease, is the most common cause of fatty liver, with a prevalence as high as 30% in many populations.¹ NAFLD may lead to fibrosis,² cirrhosis,³ liver cancer,^{4,5} liver failure requiring liver transplant,⁶ and mortality,⁷ and it is associated with type 2 diabetes, metabolic syndrome, and other cardiovascular risk factors.^{8,9} Although NAFLD represents a

Abbreviations: CI, confidence interval; CT, computed tomography; MRI, magnetic resonance imaging; MRS, magnetic resonance spectroscopy; NAFLD, nonalcoholic fatty liver disease; NIR, not reported; OR, odds ratio; QUADAS, Quality Assessment of Diagnostic Accuracy Studies; ROC, receiver operating characteristics; STARD, STAndards for the Reporting of Diagnostic accuracy studies.

From the ¹Department of Medicine, The Johns Hopkins School of Medicine, Baltimore, MD; ²Department of Medicine, Washington Hospital Center, Washington, DC; ³Department of Epidemiology, The Johns Hopkins Bloomberg School of Public Health, Baltimore, MD; ⁴Department of Radiology, The Johns Hopkins University School of Medicine, Baltimore, MD; ⁵Welch Center for Prevention, Epidemiology, and Clinical Research, The Johns Hopkins University, Baltimore, MD; ⁶Department of Cardiovascular Epidemiology and Population Genetics, National Center for Cardiovascular Research (CNIC), Madrid, Spain

This study was supported by the American Diabetes Association Mentor-Based Postdoctoral Fellowship Program (7-07-MN-08; to R.H. and M.L.), National Institute of Diabetes and Digestive and Kidney Diseases grant 1R01DK083393-01A1 (to J.M.C.), and K24-DK62222 P60 DK079637 (to F.L.B.).

*These authors contributed equally to this work.

Address reprint requests to: Ruben Hernaez, M.D., Ph.D., Department of Medicine, The Johns School of Medicine, 2024 East Monument Street, Suite 2-600, Baltimore, MD, 21287. E-mail: rhernae1@jhmi.edu; fax: 410-955-0476.

Copyright © 2011 by the American Association for the Study of Liver Diseases.

View this article online at wileyonlinelibrary.com.

DOI 10.1002/hep.24452

Potential conflict of interest: Nothing to report.

Additional Supporting Information may be found in the online version of this article.

1082

Clinically Proven Liver Steatosis Assessment with Hepatorenal Ratio



Sonographic Hepatorenal Ratio: A Noninvasive Method to Diagnose Nonalcoholic Steatosis

Valéria Ferreira de Almeida e Borges, MD,¹ Angélica L. D. Diniz, PhD,¹ Helma P. Cotrim, PhD,² Haroldo L. O. G. Rocha, MD,¹ Nestor Barbosa Andrade, MD¹

¹ Universidade Federal de Uberlândia, Uberlândia, Minas Gerais, Brazil
² Universidade Federal da Bahia, Salvador, Bahia, Brazil

Received 13 May 2011; accepted 13 August 2012

ABSTRACT: *Purpose.* To evaluate the accuracy of the sonographic hepatorenal ratio (HRR) in the diagnosis of nonalcoholic steatosis. **Keywords:** fatty liver; ultrasound; hepatorenal ratio; liver biopsy; noninvasive

Significant correlation was found between HRR and histologic steatosis ($r=0.80$). The cutoff for predicting steatosis was 1.24 (sensitivity 92.7%; specificity 92.5%).

HRR cutoff for predicting steatosis was ≥ 1.24 (sensitivity, 92.7%; specificity, 92.5%). The mean \pm SD HRRs in controls and steatosis subgroups were control 1.09 ± 0.13 , mild 1.46 ± 0.24 , moderate 1.52 ± 0.27 , severe 2.04 ± 0.3 and were significantly different from each other except between mild and moderate steatosis subgroups.

Conclusions. The HRR is a noninvasive, objective, and simple method that could be used to diagnose and grade hepatic steatosis. © 2012 Wiley Periodicals, Inc. *J Clin Ultrasound* 41:18–25, 2013; Published online in Wiley Online Library (wileyonlinelibrary.com). DOI: 10.1002/jcu.21994

Correspondence to: V. Ferreira de Almeida e Borges

© 2012 Wiley Periodicals, Inc.

18



Submit a Manuscript: <http://www.wjgnet.com/esps/>
Help Desk: <http://www.wjgnet.com/esps/helpdesk.aspx>
DOI: 10.3748/wjg.v20.i47.17985

World J Gastroenterol 2014 December 21; 20(47): 17985-17992
ISSN 1007-9327 (print) ISSN 2219-2840 (online)
© 2014 Baishideng Publishing Group Inc. All rights reserved.

CLINICAL TRIALS STUDY

Ultrasound hepatic/renal ratio and hepatic attenuation rate for quantifying liver fat content

Bo Zhang, Fang Ding, Tian Chen, Liang-Hua Xia, Juan Qian, Guo-Yi Lv

Bo Zhang, Fang Ding, Tian Chen, Liang-Hua Xia, Juan Qian, Guo-Yi Lv, Department of Echocardiography, East Hospital, Tongji University School of Medicine, Shanghai 200002, China
Guo-Yi Lv, Department of Radiology, Huazhong University of Science and Technology, Tongji Medical College, Pu Ai Hospital, Wuhan 430030, Hubei Province, China
Author contributions: Zhang B and Lv GY designed the study and wrote the manuscript; Ding F, Chen T and Xia LH performed the experiments; Qian J provided vital help in editing the manuscript.
Correspondence to: Guo-Yi Lv, MD, Deputy Chief of Physiological Department of Radiology, Huazhong University of Science and Technology, Tongji Medical College, Pu Ai Hospital, Qianzheng Road, Qianjiang District, Wuhan 430030, Hubei Province, China. laodong999@sina.com
E-mail: +86-27-68834835

Received: May 24, 2014
Accepted: September 29, 2014
Published online: December 21, 2014

liver fat content by ultrasound (quantitative ultrasound model) is: liver fat content (%) = $61.519 \times$ ultrasound hepatic/renal ratio + $167.701 \times$ hepatic echo-intensity attenuation rate - 26.736 . Spearman correlation analysis revealed that the liver fat content ratio of the quantitative ultrasound model was positively correlated with serum alanine aminotransferase, aspartate aminotransferase, and triglyceride, but negatively correlated with high density lipoprotein cholesterol. Receiver operating characteristic curve analysis revealed that the optimal point for diagnosing fatty liver was 9.15% in the quantitative ultrasound model. Furthermore, in the quantitative ultrasound model, fatty liver diagnostic sensitivity and specificity were 94.7% and 100.0%, respectively, showing that the quantitative ultrasound model was better than conventional ultrasound methods or the combined ultrasound hepatic/renal ratio and hepatic echo-intensity attenuation rate. If the ¹H-MRS liver fat content had a value < 15%, the sensitivity and specificity of the ultrasound quantitative model would

The Hepatic Renal Ratio is the strongest predictor of liver fat content.

RESULTS: Correlation analysis revealed that the ultrasound hepatic/renal ratio and hepatic echo-intensity attenuation rate were significantly correlated with ¹H-MRS liver fat content (ultrasound hepatic/renal ratio: $r = 0.952$, $P = 0.000$; hepatic echo-intensity attenuation $r = 0.850$, $P = 0.000$). The equation for predicting

Core tip: The quantitative ultrasound model is a simple, low-cost, and sensitive tool that can accurately assess hepatic fat content in clinical practice. It provides an easy and effective parameter for early diagnosis of mild hepatic steatosis and evaluation of the efficacy of

Gastrointestinal Imaging • Original Research

Hepatorenal Index as an Accurate, Simple, and Effective Tool in Screening for Steatosis

OBJECTIVE. The hepatorenal index has been reported to be a sensitive and noninvasive test to quantify steatosis, but it is cumbersome and time-consuming and requires specialized software. The aim of this study was to improve and simplify the hepatorenal index calculation and determine whether it is an effective tool for differentiating patients with steatosis from those without steatosis, thereby eliminating the need for biopsy in a large number of patients.

MATERIALS AND METHODS. One hundred one patients who had undergone ultrasound-guided percutaneous liver biopsy at our institution were selected from a patient database. Patients with renal disease, patients with liver masses, and patients whose liver and right kidney were not included on the same image were excluded. Images were acquired with high-resolution ultrasound, and the hepatorenal index was calculated using freeware based on comparison of hepatic and renal brightness.

RESULTS. Of the 101 patients, 63 had 5% or less steatosis and 38 had more than 5% steatosis.

An HRI > 1.28 has 100% sensitivity for predicting a >5% level of fat.

thirds of obese patients, and is recognized as the most common cause of chronic liver disease worldwide [1–3]. It is defined as liver fat content of more than 5% and can have no symptoms or can progress to fulminant liver failure perpetuated by a cascade of inflammatory cytokines produced by the fatty liver [4], leading to death without transplant. Steatosis is associated with multiple causes including viruses, such as the hepatitis B virus and hepatitis C virus; drugs; alcohol; idiopathic causes; obesity; and diabetes [5]. Nonalcoholic fatty liver disease (NAFLD) is becoming more common and is known to be part of a larger metabolic syndrome with potential to progress to steatohepatitis (nonalcoholic steatohepatitis [NASH]), cirrhosis, and even death [6]. In patients with diabetes, NAFLD increases the risk of hepatocellular carcinoma and has been shown to

sustained viral response rate to therapy in patients with hepatitis C [8].

Until recently, the use of ultrasound in the evaluation of steatosis has been limited because of interobserver variability of increased attenuation of sound and echogenicity in the liver and the lack of sensitivity of these findings in patients with low levels of steatosis [9]. Characterization of steatosis is typically performed using percutaneous imaging-guided biopsy, which is a costly and invasive procedure and carries the risk of adverse reactions that range from minor pain at the injection site to more serious events such as hemorrhage and death [10]. Multiple studies have sought other ways of quantifying steatosis because this important measure is useful in determining disease progression and response to therapy. Others have compared liver and kidney brightness sonographically; however, their

Clinically Proven

What the experts are saying



“The hepatorenal index is a simple, reliable, and cost-effective screening tool for identifying patients who should not undergo liver biopsy for evaluation of steatosis.”

Dr. Richard H. Marshall

Assistant Professor of Clinical Radiology
LSU School of Medicine
New Orleans, Louisiana



“The use of the hepatorenal sonographic index facilitates quantification of liver steatosis, even in small degrees, and is not affected by fibrosis or steatohepatitis. It is reproducible and operator-independent and easily can be made available and applicable in routine clinical practice.”

Dr. Muriel Webb

Department of Gastroenterology
Sourasky Medical Center
Tel Aviv, Israel

Advantages of Aixplorer®

Comparison to Fibroscan®




| Aixplorer® | FibroScan® | Aixplorer Advantage |
|--|---|--|
| <p>One Probe Solution</p> <p>XC6-1 probe works in a broad range of patients</p> | <p>Fails under common conditions of obesity and ascites¹</p> <p>Requires 3 probes with annual recalibration/replacement¹</p> | <p>Superior Workflow</p> <p>Fewer Failed Exams</p> |
| <p>60 Second Exam</p> <p>Real-time Imaging over a large 2D area</p> <p>As few as 3 acquisitions</p> <p>Qualitative and Quantitative ShearWave™ Analysis</p> <p>Steatosis Assessment with B-mode Ratio</p> <p>Fully validated connectivity options</p> | <p>5-10 Minutes per case</p> <p>Requires minimum of 10 acquisitions</p> <p>Higher technical failure rate due to blind positioning and small target area</p> <p>CAP tool adds add'l \$26k cost</p> <p>Limited connectivity</p> | <p>Higher Patient Throughput</p> <p>Lower Cost</p> |
| <p>Clinically Proven</p> <p>In HBV, HCV, NAFLD and General Liver Disease and Cirrhosis</p> | <p>Higher Technical Failure Rates</p> <p>Proven lower accuracy</p> | <p>Clinically Superior Performance and Accuracy</p> |
| <p>Higher Reimbursement</p> | <p>Low Reimbursement</p> <p>Annual Calibration Costs</p> | <p>Increased Return on Investment</p> |

¹ Elastography Assessment of Liver Fibrosis: Society of Radiologists in Ultrasound Consensus Conference Statement. Radiology: Volume 276: Number 3—September 2015

Clinically Proven SWE™ - A Solution For All Stages



SuperSonic Imagine's ShearWave™ Elastography provides solutions at ALL stages of liver disease.





Thank you!

