


PeriFlux 6000 | $tcpO_2$ made intelligent




Transcutaneous oxygen





Reveal the wound healing capacity by performing an intelligent tcpO₂

Patient born 1940
Bilateral ulcerations
Unreliable ankle pressures
Suspected ischemia



Prepare for tcpO₂
Connect electrodes and document positioning



Start tcpO₂
Follow instructions on screen

Baseline tcpO₂ value
Right foot = 1 mmHg
Left foot = 45 mmHg

What does a tcpO₂ value < 40 mmHg indicate?
Peripheral Arterial Disease (PAD)?
Capillary impairment?
Cardiopulmonary disease?
High tissue consumption of O₂ due to infection/inflammation?




Extend with a provocation
10 minutes O₂ challenge
Right foot = 1 mmHg
Left foot = 103 mmHg

Conclusions of the examination
Right foot:
Low baseline value tcpO₂. No response to O₂ challenge. If no cardiopulmonary problems are present, PAD is the most probably cause.
» Right foot critical limb ischemia with severely affected microcirculation
Left foot:
Borderline tcpO₂ baseline value. Good response to O₂ challenge.
» Beneficial conditions for spontaneous healing

tcpO₂ to predict benefit from hyperbaric oxygen treatment

Hyperbaric oxygen therapy (HBOT) is a medical treatment in which the patient breathes 100% oxygen inside a pressurized chamber. The aim is to increase the general oxygen level in the tissue. It is commonly used in conditions where partial or total tissue hypoxia is present, such as diabetic foot lesions or arterial ulcers. Not all patients will benefit from HBOT. Therefore, it is often recommended to use tcpO₂ to qualify patients before initiating the therapy. Furthermore, tcpO₂ can monitor the tissue oxygenation throughout the treatment. Perimed offers a remote panel system to facilitate measurements inside a multi-place chamber. In addition, all equipment is compatible with extension cables and penetrators used in mono-chambers.



Microvascular assessment
OPERATOR: Anna Bauer DATE: 2012-09-14 TIME: 09:29 HOSPITAL LOGO Vascular Lab Department

PATIENT INFORMATION
NAME: John Smith
ID: 1234
DATE OF BIRTH: 1940-05-20
GENDER: Male

SUMMARY / CONCLUSION
Right foot CLI with severely affected microcirculation. Left foot beneficial conditions for healing.

ANAMNESIS
DIABETES: Yes
INSULINE DEPENDANT: No
DIABETES SINCE: 1998
NEUROPATHY: YES
HYPERTENSION: Yes
MEDICALLY REGULATED: No
BLOOD PRESSURE: 160/70
SMOKER: Yes
PACK PER YEAR: 1500

LEFT
WOUNDS: Yes
WALKING PAIN: Yes
REST PAIN: Yes
PALMABLE PULSES: No

RIGHT
WOUNDS: Yes
WALKING PAIN: Yes
REST PAIN: Yes
PALMABLE PULSES: Yes

DETAILS:
Wound 6 months

TCOM RESULTS

NAME	BL (mmHg)	O ₂ (Pa)
1. Site 1	1	1 (0%)
2. Site 2	2	3 (+50%)
3. Site 3	1	1 (+0%)
4. Site 4	38	82 (+110%)
5. Site 5	45	103 (+129%)
6. Site 6	47	105 (+123%)
7. Site 7	43	102 (+137%)
8. Site 8	49	105 (+114%)

Recorded using PeriFlux 6000

PeriFlux 6000 - tcpO₂ / tcpCO₂ monitoring

Why measure tcpO₂?

Peripheral Arterial Disease (PAD) will influence the outcome of wound healing. To prevent limb loss, the degree of ischemia has to be assessed at an early stage. Transcutaneous oxygen (tcpO₂) has proven to be a useful tool for this purpose as it reflects the metabolic state of the lower limb. tcpO₂ is particularly important for the assessment of wounds and prediction of amputation levels in patients with critical limb ischemia and/or diabetes, as these patients commonly also have impaired microvascular function and falsely elevated ABI values.

What is tcpO₂ and tcpCO₂?

Transcutaneous oxygen, tcpO₂ or TCOM, is a local non-invasive measurement reflecting the amount of O₂ that has diffused from the capillaries, through the epidermis, to a Clark-type electrode at the measuring site. It provides instant, continuous information about the body's ability to deliver oxygen to the tissue. Any impairment in the ability to deliver oxygen to the tissue will be revealed immediately since the skin is ranked very low in the body's system of oxygenation priority. tcpO₂ measurements usually require at least 2-3 sites, preferably 4 or more, to provide a good picture of the oxygenation status of the skin.



Transcutaneous carbon dioxide (tcpCO₂) is the local carbon dioxide tension in the skin. The carbon dioxide level is affected by both the local metabolic processes and the blood perfusion ability to remove carbon dioxide.

Reference values tcpO₂

50-70 mmHg Normal
< 40 mmHg Impaired wound healing
< 30 mmHg Critical Limb Ischemia

Interpreting results

As a tcpO₂ value can be influenced by many factors, extending the baseline with provocations, or including a reference electrode, may prove valuable for the clinical decision making.

Oxygen challenge (tcpO₂ measurement during 100 % oxygen inhalation) will distinguish low values due to a barrier to oxygen diffusion (edema and/or inflammation) from macrovascular disease (PAD).

Leg elevation for a duration of 5-15 minutes may be used to confirm macrovascular disease. Other methods to confirm macrovascular disease include toe and ankle pressure.

Reference electrode or oxygen saturation (pulse oximeter) will rule out arterial hypoxemia (due to pulmonary disease, for example).



PeriFlux 6000 is the latest generation transcutaneous O₂ and CO₂ equipment from Perimed



PeriFlux 6000 Specifications

Start-up time: Maximum 60 seconds
 Automatic calibration: In air (tcpO₂) / with TC 600 (tcpCO₂), 8 electrodes simultaneously
 Memory storage capacity: 2 GB
 Alarm: Visual and audible
 Dimensions: W=28 cm, H=22 cm, D=25 cm
 Weight: 4.9 kg (equipped with 8 PF 6040 units)
 Display: Touchscreen: 8.4" color TFT-LCD, Resolution: 800x600 px
 Power consumption: 100 to 240 VAC, 50 or 60 Hz, 65 VA
 Operating conditions: Temp.: +15 to +35 °C at 10 to 85 % RH, Environmental pressure: 70 to 110 kPa / 700 to 1100 mbar
 External connections: 2 USB hosts (for connecting printer, camera, keyboard, pointer device, etc.), 1 USB device (for connecting PC)
 Humidity sensor: Range: 10 to 85 % RH, Accuracy: ± 4 % RH

PF 6040 tcpO₂/tcpCO₂ Unit

One electrode per unit
 Measured parameters: tcpO₂, tcpCO₂
 Measurement ranges: tcpO₂ = 0 to 1999 mmHg (0–267 kPa), tcpCO₂ = 5 to 200 mmHg (0.67–26.7 kPa)
 Accuracy: tcpO₂ < ±5 mmHg from 0 to 20.9 % O₂ and < ±10 % of reading from 21% to full scale
 tcpCO₂ ±5 mmHg over measurement range (5 to 100 mmHg)
 Temperature settings: Range: 37 to 45 °C, set in steps of 0.5 °C, Accuracy: 0.5 °C
 Built-in barometer: Range: 225 to 825 mmHg, Accuracy: ± 3.0 mmHg
 Classification type: BF (body floating)

Electrodes:

E5250: pO₂ sensor
 E5280: Combined pO₂ / pCO₂ sensor

Compliance:

HIPAA compliant
 MDD 93/42/EEC, WEEE 2002/96/EG, ROHS 2002/95/EG, EN60601-1:2006 (Third edition), EN60601-1-2:2007, EN60601-1-6:2010, ASTM D4169:2009, EN ISO10993-1:2009, EN62304:2006, 21 CFR 800-1299:2008, ANSI/AAMI ES60601-1:2005, CMDR, 2010, CAN/CSA-C22.2 No. 60601-1-08, IEC60601-2-23:2011, EN60601-1-8:2007 (Second edition), NFPA 99:2012, GB 18455-2001, SJ/T 11363-2006, SJ/T 11364-2006, EN 980:2008, ISO15223-1:2007 (First edition), EN62366:2008, EN 1041:2008, MEDDEV. 2.7.1 Rev.3, EN ISO 14971:2012

Accessories and Consumables:

Fixation rings: TC 550 Fixation Rings for tcpO₂ / tcpCO₂, TC 555 Fixation Rings Extra Strength Adhesive for tcpO₂ / tcpCO₂
 Contact liquid (20 ml): TC 560 Contact Liquid
 Membraning kit: D826 Membraning Kit tcpO₂, D280 Membraning Kit tcpCO₂
 Calibration unit for CO₂: TC 600 Calibration Unit
 Calibration gas (CO₂): TC 510 Calibration Gas
 Remote panel: PF 5840 TC Remote Panel
 Cables for remote panels: PF 5841 Extension Cable 3 m, PF 5842 Extension Cable 6 m
 Color coded labels: PF 6103 Color Coded Labels
 Camera: PF 6113 Camera
 System carts
 Foot pedal
 Demand valve EASE II 03 3M SS/DIN 120 and range of different sized masks
 Medical isolation transformer, Network isolator

Due to Perimed's commitment to continuously improve our products, all specifications are subject to change without notice. The 510(k) approval for the PeriFlux 6000 does not yet cover the modules PF 6010 and PF 6050.

Standard PeriFlux 6000 configurations:

	tcpO ₂	PRESSURE Standard	COMBINED Standard	PRESSURE Premium	COMBINED Premium
Toe pressure, ABI and PVR	-	●	●	●	●
tcpO ₂	1 - 8	-	2	-	3
Treadmill	-	○	○	●	●
Segmental pressure	-	○	○	●	●


○ Available ● Included - Not applicable

For more information please contact Perimed AB

Perimed AB, Datavägen 9A, SE-175 43 Järfälla-Stockholm, Sweden | Tel: +46-8-580 119 90 Fax: +46-8-580 100 28
 E-mail: mail@perimed-instruments.com | Website: www.perimed-instruments.com
 Part. No. 44-00284-03



www.perimed-instruments.com



1 - 2 - 3



Compact and elegant solution

The PeriFlux 6000 is operated using a touch screen interface. It can be equipped with up to 8 channels of tcpO₂ allowing for accurate mapping of the extremity. It is small, portable and can be mounted on an arm or stand.

Step-by-step-instructions

The user is guided throughout the procedure by simple instructions displayed on the screen. Different tests may be implemented, including exam room and in-chamber measurements.

Automatic report generator

All test results, including the site positioning photo, are displayed in an automatically generated report that may be printed or exported as a PDF file. The report template can be customized according to the requirements of the user.

HIPAA compliant

For patient security the PeriFlux 6000 is HIPAA compliant.

Billing and Reimbursement

Use CPT codes 93922 and 93923 for billing and reimbursement of tcpO₂ measurements.

Connection to PC

Data from the PeriFlux 6000 can be transferred to a PC. The PeriFlux Configuration Software (PCS) also makes it possible to review generated data on a PC.

References

1. European Society for Vascular Surgery, CLI Guideline Committee Guidelines for Critical Limb Ischaemia and Diabetic Foot, 2011
2. Transcutaneous Oximetry in Clinical Practice: Consensus statements from an expert panel based on evidence. Fife, Smart, Sheffield, Hopf, Hawkins and Clarke. UHM 2009, Vol. 36, No. 1.