

VitaLogik™ 6000 Joins Mennen Medical Monitor Solution

The VitaLogik™ 6000 is a member of the VitaLogik™ family together with Envoy™, VitaLogik™ 5000 series, VitaLogik™ 4000 series and the Enguardview™ providing a complete range of patient monitoring solutions.

Making the Connection

The LAN connectivity provides bed to bed communication, alarm watch, network laser printer and connectivity to Enguardview™ Workstation, Enguard™ Remote Monitor or Ensemble™ Central Nurse Station.

Wireless Networking - For Patient Care Anywhere

The wireless network option extends Central Station monitoring and automated vital sign collection beyond the bedside.

MenneNet

The MenneNet™ (Mennen Network) comprehensive solution, addresses the increasing demands for integration and connectivity of all Mennen monitoring devices using the LAN, WLAN and WEB connectivity in the hospital.

Mennen Medical Monitor Web Portal

Mennen Medical's Monitor Web Portal is an advanced proprietary software application on the Mennen Server that provides clinicians remote access to patient data from any location using the hospital's Web intranet/internet connection. The Web Portal software application displays real time monitored patients and previously "Saved & Discharged" patients in remote view.

Envoy Explorer

Envoy Explorer Software for connectivity between a PC to the MenneNet™ enables viewing and storing of vital signs data from any Envoy™/VitaLogik™ Monitor on the network.

CE 0473 FDA Pending



www.mennenmedical.com

936-OPT-113 Rev A 12/2009

Headquarters:
Mennen Medical Ltd.
4 Ha-Yarden Street
Yavne, P.O.Box 102
Rehovot 6100, Israel
Phone: +972-8-9323333
Fax: +972-8-9328510

United States:
Mennen Medical Corp.
950 Industrial Blvd.
Southampton,
PA 18966-4070, US
Phone: (215) 259-1020
Fax: (215) 675-6212

United Kingdom:
Mennen Medical UK.
Unit 18 Avant Business Center
21 Denbigh Road
Milton Keynes MK1 1DT, UK
Phone: +44 1908 648070
Fax: +44 1908 646030

2009 Mennen Medical Ltd. All rights reserved. Specification subject to be changed without notice.

© 2009 Mennen Medical Ltd.



VitaLogik™ 6000

Configured Patient Monitor



We Follow Every Beat of Your Heart

The VitaLogik™ 6000 is a pre-configured multi-parameter Patient Monitor with a 15" flat screen display and a battery backup. The Monitor is suitable for use in a broad range of patient applications for adults, pediatrics and neonates with the following functions: display of waveforms and vital signs, alarm with status messages and event strips printing. The VitaLogik™ 6000 has an extensive storage capability of Full Disclosure, charts, trends and event strips. Data of patients discharged from the monitor is saved, for up to 10 patients per monitor.

The Vital Signs Monitor Includes:


- ECG/Heart Rate
- ST Segment Analysis
- Full Range of Arrhythmias
- Respiration
- Pulse Oximetry (Nellcor®/Masimo®)
- Non-Invasive Blood Pressure
- 2 Temperatures

Optional Additional Features:


- 2/4 Invasive Blood Pressure Channel and Cardiac Output
- Microstream® EtCO₂
- BIS Module with 2 EEG waveforms
- Anesthesia Multi Gas Module
- 3 Channel Recorder/Printer
- Touch Screen Operation
- 2nd Li-Ion Battery
- Extended Memory
- Wireless Networking
- Interface to Ventilators / Anesthesia Machine
- Interface to Vigilance™ SvO₂/CCO Monitor



Very Big Numbers (VBN)




External Device Interface



Full Disclosure Storage



12 Lead Diagnostic ECG Display



Connectivity via CITRIX



Connectivity to HIS/LIS

